

The Use of Absolute Cerebral Oximetry in Cardiovascular Surgery

Gregory W Fischer, M.D., Cardiac Anesthesiologist, Department of Anesthesiology
Konstadinos A Plestis, M.D, Cardiothoracic Surgeon, Department of Cardiothoracic Surgery,
Mount Sinai Medical Center, New York, NY, USA

Introduction

Cerebral oximetry, based on near infrared spectroscopy (NIRS) technology, provides information on the availability of oxygen in brain tissue at risk during numerous pathological conditions.¹ Cerebral oximetry measures local concentrations of hemoglobin (oxy- and deoxy-), and regional cerebral tissue oxygen saturation (SctO₂) at the microvascular level (arterioles, venules, and capillaries only).^{2,3,4} As a result, cerebral oximetry SctO₂ is a mixed oxygen saturation parameter which has a value between arterial (SaO₂) and jugular venous oxygen saturation (SjvO₂) under normal physiological conditions, therefore SaO₂ > SctO₂ > SjvO₂. Complementary to the arterial oxygen saturation (SaO₂) measured by pulse oximetry, SctO₂ reflects regional cerebral metabolism and the balance of local cerebral oxygen supply/demand. The advantages of cerebral oximetry are: 1) It provides SctO₂ values continuously and non-invasively at the bedside;⁵ 2) SctO₂ is a sensitive index of cerebral hypoxia and/or cerebral ischemia^{6,7} which is one of the main causes of brain injury in clinical settings^{8,9}.

The FORE-SIGHT™ Cerebral Oximeter (CAS Medical Systems) is significantly different from cerebral oximeters currently on the market. The FORE-SIGHT monitor was developed with the support of a series of Small Business Innovation Research Grants from the National Institute of Neurological Disorders and Stroke (NINDS) of the National Institute of Health (NIH).¹⁰ It is the only absolute cerebral oximeter cleared by the FDA¹¹ based on accuracy. The FORE-SIGHT Cerebral Oximeter, with its ability to provide absolute measurement makes it possible to establish threshold values for SctO₂ that can be used to guide clinical interventions.

FORE-SIGHT Cerebral Oximeter determined cerebral tissue oxygen saturation, SctO₂, is defined as the ratio of concentrations of HbO₂ and Hb + HbO₂ in the brain tissue, thus $SctO_2 = 100\% * HbO_2 / (Hb + HbO_2)$. The value of SctO₂ reflects a proportional mix of arterial and venous blood that can be calibrated from arterial and internal jugular venous blood¹². It is estimated that the NIRS cerebral oximeter interrogated brain tissue microvasculature is about 70% venous and 30% arterial during most physiological conditions in humans based on Positron Emission Tomography (PET) studies¹³. In validation studies, FORE-SIGHT cerebral oximeter determined SctO₂ showed a strong correlation with the reference SctO₂ over the spectrum of pulse oximeter determined arterial oxygen saturation SpO₂ values between 70 and 100% from 18 subjects.^{14,15} The bias and precision (1 standard deviation) for the FORE-SIGHT Cerebral Oximeter SctO₂ compared to reference SctO₂ derived from co-oximetry of arterial and jugular bulb blood was 0.18 ± 3.7 (1SD). The FORE-SIGHT Cerebral Oximeter, with its ability to provide absolute measurement of cerebral tissue oxygen saturation SctO₂, overcomes the limitations of previous cerebral oximeters.¹⁶ Particularly, clinically relevant SctO₂ threshold values can be established with the FORE-SIGHT Cerebral Oximeter for physicians to provide tailored patient management¹⁷. It is known that SjvO₂ has a normal lower limit at ~45% and the upper limit at 70%.^{18,19} The FORE-SIGHT SctO₂ is about 10% higher than SjvO₂ consistently over a wide range of oxygen saturation values. Therefore, the absolute FORE-SIGHT Cerebral Oximeter lower safe SctO₂ threshold is about 55%.

The Need for Bedside Cerebral Oximetry

Monitoring brain oxygenation is critical in providing information used to guide patient management in many clinical situations.^{20,21} Currently, brain oxygenation can be measured invasively by jugular bulb oximetry SjvO₂ or brain tissue pO₂ sensor.^{22,23} Benefits of SjvO₂ monitoring include: a) improved outcome in physiologic management of head injury patients (370,000 cases/year, in the US);^{24,25,26} b) detection of critical

events when brain oxygenation could be compromised during cardiac surgery (800,000 cases/year worldwide)²⁷, and neurosurgery.²⁸ In some institutions, S_{ijv}O₂ monitoring is routinely applied in surgery of the aorta²⁹, and in neuro intensive care units.³⁰ Therefore, a bedside cerebral oximeter that can provide non-invasive measurement of cerebral oxygenation is highly desirable.

Despite decline in overall mortality after coronary artery bypass grafting (CABG) and valvular surgery with cardiopulmonary bypass (CPB), the rates of cognitive dysfunction have not improved.³¹ In some reports, most patients had subtle signs of impaired cognitive performance, with incidences ranging from 60 to 80%.³² There are two different forms of brain injury that may occur after CABG or thoracic aorta surgery: neurological dysfunction (ND) and neurocognitive dysfunction (NCD). ND is defined as clinically evident focal or global neurological injury resulting in stroke, hypoxic encephalopathy, transient ischemic attack, or stupor. NCD is defined as postoperative confusion, agitation, delirium, prolonged obtundation, or transient Parkinsonism. NCD occurs more frequently, affecting 40% to 80% of CABG/aortic surgery patients, depending on the method of detection.³³ While it is easy to diagnose post-operative ND, NCD is more subtle, and needs to be evaluated with a full battery of neurocognitive testing administered by trained professionals. Recent reports based on accurate neurocognitive testing before and after surgery have suggested that NCD can no longer be considered a benign self-limiting condition, but rather a long-lasting neurocognitive insult capable of reducing quality of life by impairing memory and fine motor function.³⁴

The etiology of brain injury following cardiac surgery is still not completely understood and somewhat controversial. Some of the possible mechanisms include diffuse microembolization, cerebral hypoperfusion, and metabolic factors; the incidence of injury seems to be higher when CPB duration exceeds 70 minutes, and when there is rapid rewarming, particularly in the older surgical population.³⁵ Regardless of the immediate cause, such persistent cognitive dysfunction likely results from brain ischemia during surgery, which may be a result of focal arterial embolism,³⁶ global hypo-perfusion of the brain,³⁷ or an interaction of the two.³⁸

Strategies for preventing arterial embolism and brain hypo-perfusion differ. To avoid arterial embolism during CPB, arterial line filters, intraoperative imaging, and careful manipulation by the surgeon are essential. On the other hand, diffuse hypo-perfusion of the brain can be avoided only by very careful planning: if something goes wrong, the only hope is early detection and immediate restoration of adequate perfusion before irreversible brain damage develops. For this purpose, sensitive, real-time monitoring of brain ischemia during such surgical procedures is needed.³⁹ At the present time cerebral oximetry is the only feasible technology that monitors cerebral hypoxia and/or cerebral ischemia noninvasively and continuously.

Cerebral oximeters provide information that other bedside brain monitors, such as electroencephalography (EEG) and transcranial Doppler (TCD), cannot offer. Other modalities, such as positron emission tomography (PET), perfusion CT, and magnetic resonance imaging (MRI) can provide detailed “snap shot” information about cerebral oxygenation, but cannot be used at the bedside. Since cerebral oximetry can provide an immediate indication of cerebral blood flow changes, oxygenation changes, it could find a wide range of applications in operating rooms (OR), in recovery rooms, as well as in intensive care units (ICU). In all of these situations, cerebral oximetry can be used to monitor the safety and efficacy of treatment interventions.

Why Absolute Cerebral Oximetry

Current trend only cerebral oximetry may be inadequate:

Previous cerebral oximeters on the market measure cerebral oxygenation as a trend only.⁴⁰ Some studies suggest that operative technique can be modified based on application of trend-only cerebral oximetry. Since these monitors measure trends only, a baseline first has to be established and cerebral oximetry values need to be maintained at or near preoperative baseline.^{41,42} Another approach is to keep the cerebral oxygen saturation at levels within 20-25% of the anesthesia pre-induction value.^{43,44} However, studies have shown that seventy five percent of patients undergoing coronary bypass have a significant impairment in baseline regional cerebral perfusion (rCP).⁴⁵ Other studies also demonstrated that patients undergoing cardiovascular

surgery have a high prevalence of cerebral vascular disease in varying degrees.⁴⁶ In addition, abnormal preoperative rCP was found to be a strong indicator for post surgical decline in neuropsychologic testing.⁴⁷ These findings suggest that it is difficult to define a “normal” pre-induction baseline value for the trend only cerebral oximetry. A percentage drop based on the unreliable baseline value is more questionable to serve as a threshold for clinical intervention. In fact, studies have confirmed that while a trend only cerebral oximeter can detect adverse brain oxygenation by measuring the change in SctO₂ from a baseline value, it cannot provide accurate and reliable normal and threshold values of cerebral tissue oxygen saturation.^{48,49,50,51,52,53,54}

Absolute cerebral oximetry is essential for tailored patient management:

There is an increasing amount of evidence that has demonstrated the need for a tailored patient management protocol as current approaches for managing flow, arterial blood pressure, and pH during cardiac surgery are based on studies that included few elderly or high-risk patients and predated many other contemporary practices⁵⁵. For example, watershed-distribution stroke happens more frequently in patients undergoing cardiac surgery than in general stroke population (over 40% versus 2-5%, respectively).^{56,57,58} Gottesman et. al., reported that mechanism of watershed stroke after cardiac surgery may include an intraoperative drop in blood pressure from a patient’s baseline.⁵⁹ This suggests that following the standard protocol to maintain an optimal range for blood pressure during cardiac surgery is insufficient for some patients. We believe that what is needed, is an online monitor to evaluate the effect of blood pressure level as well as blood pressure change on the brain. Our own studies suggest that the FORE-SIGHT Absolute Cerebral Oximeter could be used for this tailored patient management approach. Readings of cerebral tissue oxygen saturation SctO₂ indicated that maintaining mean arterial blood pressure at 50-60 mmHg during hypothermic CPB is tolerated by most patients, but this level seems to be inadequate for certain patients.

References:

- ¹ Ferrari M, Mottola L, Quaresima Principles, techniques, and limitations of near infrared spectroscopy. *V.Can J Appl Physiol.* 2004 Aug;29(4):463-87.
- ² Benni PB, Chen B, Dykes FD, Wagoner SF, Heard M, Tanner AJ, Young TL, Rais-Bahrami K, Rivera O, Short BL, “Validation of the CAS Neonatal NIRS System by Monitoring VV-ECMO Patients: Preliminary Results”, *Experimental Medicine and Biology*, 2005;Volume 566.
- ³ Kurth CD, et. al., Cerebral oxygen saturation before congenital heart surgery. *Ann Thorac Surg.* 2001 Jul;72(1):187-92.
- ⁴ Watzman, HM, Kurth C.D, Montenegro LM, et al. Arterial and venous contributions to near infrared cerebral oximetry. *Anesthesiology* 2000;93:947–53.
- ⁵ Madsen, P.L. and N.H. Secher, Near-infrared oximetry of the brain. *Prog Neurobiol*, 1999. 58(6): p. 541-60.
- ⁶ Casati A, Fanelli G, Pietropaoli P, Proietti R, Tufano R, Danelli G, Fierro G, De Cosmo G, Servillo G. Continuous monitoring of cerebral oxygen saturation in elderly patients undergoing major abdominal surgery minimizes brain exposure to potential hypoxia, *Anesth Analg.* 2005 Sep;101(3):740-7
- ⁷ Yao FS, Tseng CC, Ho CY, Levin SK, Illner P. Cerebral oxygen desaturation is associated with early postoperative neuropsychological dysfunction in patients undergoing cardiac surgery. *J Cardiothorac Vasc Anesth.* 2004 Oct;18(5):552-8.
- ⁸ Orihashi K, Sueda T, Okada K, Imai K. Eur Near-infrared spectroscopy for monitoring cerebral ischemia during selective cerebral perfusion. *J Cardiothorac Surg.* 2004 Nov;26(5):907-11.
- ⁹ Hogue CW Jr, Palin CA, Arrowsmith JE. Cardiopulmonary bypass management and neurologic outcomes: an evidence-based appraisal of current practices. *Anesth Analg.* 2006 Jul;103(1):21-37.
- ¹⁰ <http://crisp.cit.nih.gov/>
- ¹¹ www.casmed.com
- ¹² McCormick PW, et. al. Noninvasive cerebral optical spectroscopy for monitoring cerebral oxygen delivery and hemodynamics. *Crit Care Med* 1991 Jan;19(1):89-97
- ¹³ Ito H, Kanno I, Fukuda H; Human cerebral circulation: positron emission tomography studies. *Ann Nucl Med.* 2005 Apr;19(2):65-74. Review.
- ¹⁴ MacLeod DB, Ikeda K, Keifer JC, Moretti E, and Ames W, Validation of the CAS Adult Cerebral Oximeter during Hypoxia in Healthy Volunteers, *Anesth Analg* 2006; 102:S162
- ¹⁵ MacLeod DB, Ikeda K, Moretti E, Keifer JC, and Grocott H, Using the CAS Cerebral Oximeter to Estimate Cerebral Venous Oxygen Saturation Presented at the ASA 2005

-
- ¹⁶ Taillefer MC, Denault AY. Cerebral near-infrared spectroscopy in adult heart surgery: systematic review of its clinical efficacy. *Can J Anaesth*. 2005 Jan;52(1):79-87.
- ¹⁷ Fischer GW, Reich D, Plestis KA, Griep RB, Results Utilizing Absolute Cerebral Oximetry Monitoring Suggests the Need for Tailored Patient Management during Cardiac Surgery, Presented at the Outcomes 2006: "The Key West Meeting".
- ¹⁸ Macmillan CSA and Andrews PJD, Cerebrovenous oxygen saturation monitoring: practical considerations and clinical relevance. *Intensive Care Med* 2000; 26:1028-36
- ¹⁹ Chierigato A, et. al., Normal jugular bulb oxygen saturation. *J Neurol Neurosurg Psychiatry* 2003; 74:784-786.
- ²⁰ Feldman Z, et. al., Monitoring of cerebral hemodynamics with jugular bulb catheters. *Crit Car Clin* 1997; 13:51-77.
- ²¹ De Deyne C, et. al., Jugular bulb oximetry: review on a cerebral monitoring technique. *Acta Anaesthesiol Belg* 1998; 49:21-31.
- ²² al-Rawi PG, et. al., Multiparameter brain tissue monitoring--correlation between parameters and identification of CPP thresholds. *Zentralbl Neurochir* 2000;61(2):74-9
- ²³ Jodicke A, et. al., Monitoring of brain tissue oxygenation during aneurysm surgery: prediction of procedure-related ischemic events. *J Neurosurg* 2003 Mar;98(3):515-23
- ²⁴ Cruz J. The first decade of continuous monitoring of jugular bulb oxyhemoglobin saturation: Management strategies and clinical outcome. *Cri Care Med* 1998; 26:344-51.
- ²⁵ Gopinath Sp. et. al., Jugular venous desaturation and outcome after head injury. *J neurol Neurosurg Psychiatr* 1994; 57:717-23.
- ²⁶ Cormio M, et.al., Elevated jugular venous oxygen saturation after severe head injury. *J Neurosurg* 1999; 90: 9-10
- ²⁷ Croughwell ND, et. al., Jugular bulb saturation and cognitive dysfunction after cardiopulmonary bypass. *Ann Thorac Surg* 1994; 58:1072-8.
- ²⁸ Moss E, et. al., Effects of changes in mean arterial pressure on S_jO₂ during cerebral aneurysm surgery. *Br J Anaesth* 1995; 75:527-30.
- ²⁹ Reich DL, et. al. Using jugular bulb oxyhemoglobin saturation to guide onset of deep hypothermic circulatory arrest does not affect post-operative neuropsychological function. *Eur J Cardiothorac Surg*. 2004 Mar;25(3):401-6; discussion 406-8.
- ³⁰ Himmelseher S, et. al., Intraoperative monitoring in neuranesthesia: a national comparisons between two surveys in Germany in 1991 and 1997. *Anesth Analg* 2001; 92: 166-71
- ³¹ Roach GW, et al. Adverse cerebral outcomes after coronary bypass surgery. *N Engl J Med* 1996;335:1857-63.
- ³² Newman MF, et al. Longitudinal assessment of neurocognitive function after coronary artery bypass surgery. *N Engl J Med* 2001;344:395-402.
- ³³ Richard FM, Bill IW: Normothermic versus hypothermic cardiopulmonary bypass: Central nervous system outcome. *J Cardiothorac Vasc Anesth* 10:45-53, 1996
- ³⁴ Newman MF, et al. Report of the substudy assessing the impact of neurocognitive function on quality of life 5 years after cardiac surgery. *Stroke*. 2001; 32: 2874-2881.
- ³⁵ Newman MF, et al. Longitudinal assessment of neurocognitive function after coronary artery bypass surgery. *N Engl J Med* 2001;344:395-402.
- ³⁶ Ergin MA, et. al., Hypothermic circulatory arrest in operations on the thoracic aorta. Determinants of operative mortality and neurologic outcome. *J Thorac Cardiovasc Surg* 1994; 107: 788
- ³⁷ Plestis KA, Gold JP. Importance of blood pressure regulation in maintaining adequate tissue perfusion during cardiopulmonary bypass. *Semin Thorac Cardiovasc Surg*. 2001 Apr;13(2):170-5.
- ³⁸ Christian Hagl et. al., Hypothermic circulatory arrest during ascending and aortic arch surgery: the theoretical impact of different cerebral perfusion techniques and other methods of cerebral protection *Eur J Cardiothorac Surg* 2003;24:371-378
- ³⁹ Murkin JM. Hemodynamic changes during cardiac manipulation in off-CPB surgery: relevance in brain perfusion. *Heart Surg Forum*. 2002;5(3):221-4.
- ⁴⁰ FDA 510 (K) - K001842
- ⁴¹ Goldman S, et. al., Interventions based on cerebral oximetry reduce the incidence of prolonged ventilation and hospital stay in cardiac surgery patients. *Outcomes 2004:The Key West Meeting, Florida, May 19-24, 2004*
- ⁴² Goldman S, et. al., Optimizing intraoperative cerebral oxygen delivery using noninvasive cerebral oximetry decreases the incidence of stroke for cardiac surgical patients. *Heart Surg Forum*. 2004;7(5):E376-81.
- ⁴³ Murkin JM. Perioperative detection of brain oxygenation and clinical outcomes in cardiac surgery. *Sem Cardioth Vasc Anesth* 2004; 8(1): 13-14.
- ⁴⁴ Murkin JM; et. al., Brain oxygenation in diabetic patients during coronary surgery: A randomized prospective blinded study. *Anesth Analg* 2005;100:SCA101.

-
- ⁴⁵ Moraca R, et. al., Impaired baseline regional cerebral perfusion in patients referred for coronary artery bypass. *J Thorac Cardiovasc Surg.* 2006 Mar;131(3):540-6
- ⁴⁶ Nakamura Y. et. al., The prevalence and severity of cerebralvascular disease in patients undergoing cardiovascular surgery. *Ann Thorac Cardiovasc Surg* 2004; 10:81-4
- ⁴⁷ Hall RA, et. al., Brain SPECT imaging and neuropsychological testing in coronary artery bypass patients: single photon emission computed tomography. *Ann Thorac Surg.* 1999 Dec;68(6):2082-8.
- ⁴⁸ Reents W, et. al., Cerebral oxygen saturation assessed by near-infrared spectroscopy during coronary artery bypass grafting and early postoperative cognitive function. *Ann Thorac Surg.* 2002 Jul;74(1):109-14.
- ⁴⁹ Oki A, et.al., Simultaneous monitoring of somatosensory evoked potentials and regional cerebral oxygen saturation combined with serial measurement of plasma levels of cerebral specific proteins for the early diagnosis of postoperative brain damage in cardiovascular surgery. *J Artif Organs.* 2004;7(1):13-8.
- ⁵⁰ Yoshitani K et. al., Comparison of changes in jugular venous bulb oxygen saturation and cerebral oxygen saturation during variations of haemoglobin concentration under propofol and sevoflurane anaesthesia. *Br J Anaesth.* 2005 Mar;94(3):341-6. Epub 2004 Dec 10. Comment in: *Br J Anaesth.* 2005 Jun;94(6):863; author reply 863-4.
- ⁵¹ Schwarz G, Cerebral oximetry in dead subjects. *J Neurosurg Anesthesiol.* 1996 Jul;8(3):189-93. Comment in: *J Neurosurg Anesthesiol.* 1997 Apr;9(2):194-5. *J Neurosurg Anesthesiol.* 1997 Jan;9(1):76-7.
- ⁵² Ulrich Beese et. al., Comparison of Near-Infrared Spectroscopy and Somatosensory Evoked Potentials for the Detection of Cerebral Ischemia During Carotid Endarterectomy. *Stroke.* 1998;29:2032-2037
- ⁵³ K. Büchner et. al., Near-infrared spectroscopy - not useful to monitor cerebral oxygenation after severe brain injury. *Zentralbl Neurochir* 2000; Vol. 61: 69-73
- ⁵⁴ Mille T, Near infrared spectroscopy monitoring during carotid endarterectomy: which threshold value is critical? *Eur J Vasc Endovasc Surg.* 2004 Jun;27(6):646-50.
- ⁵⁵ Hogue CW Jr, et. al., Cardiopulmonary bypass management and neurologic outcomes: an evidence-based appraisal of current practices. *Anesth Analg.* 2006 Jul;103(1):21-37.
- ⁵⁶ Yamamoto Y, Georgiadis AL, Chang H-M, Caplan LR. Posterior cerebral artery territory infarcts in the New England medical center posterior circulation registry. *Arch Neurol.* 1999;56:824–832.
- ⁵⁷ Rankin JM, et. al., Mechanism of stroke complicating cardiopulmonary bypass surgery. *Aust N Z J Med.* 1994;24:154 – 160.
- ⁵⁸ Paciaroni M, et. al., Neurovascular territory involved in different etiological subtypes of ischemic stroke in the Perugia stroke registry. *Eur J Neurol.* 2003;10:361–365.
- ⁵⁹ Gottesman RF. et. al., Watershed Stroke After Cardiac Surgery - Diagnosis, Etiology, and Outcome, *Stroke.* 2006; 37:2306-2311