


# Tuff-Cuff® & Safe-Cuff®



Your requirements for cuffs vary by care area and patient need. CASMED® cuffs may be ordered with connectors that make them compatible with most brands of automated blood pressure monitors and wall or portable manometers.

Durable, leak-proof, comfortable, accurate, and economical

Three styles are available:

**Safe-Cuff**

- Soft fabric prevents perspiration under cuff

**Safe-Cuff Vinyl**

- Soft vinyl is comfortable yet can be wiped off

**Tuff-Cuff Reusable**

- PVC/DEHP free
- Anti-microbial protection built into the fabric
- Fluid barrier

**Features**

- Meets AHA/ AAMI Standards
- Color coded
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Latex free
- Fabric fold at elbow is more comfortable than a weld
- CASMED adapters allow for compatibility with monitors worldwide



**Fittings**



DM



HP



ML



FL



MQ



MQML



MQMF (MQ + MQML)

Description	Single tube	Dual tube	Limb circumference
Tuff-Cuff Infant	TCS0813--	TCD0813--	08 – 13 cm
Tuff-Cuff Child	TCS1219--	TCD1219--	12 – 19 cm
Tuff-Cuff Small Adult	TCS1725--	TCD1725--	17 – 25 cm
Tuff-Cuff Adult	TCS2333--	TCD2333--	23 – 33 cm
Tuff-Cuff Adult Long	TCS2333L--	TCD2333L--	23 – 33 cm
Tuff-Cuff Large Adult	TCS3140--	TCD3140--	31 – 40 cm
Tuff-Cuff Large Adult Long	TCS3140L--	TCD3140L--	31 – 40 cm
Tuff-Cuff Thigh	TCS3850--	TCD3850--	38 – 50 cm
Safe-Cuff Infant	SCS0813--	SCD0813--	08 – 13 cm
Safe-Cuff Child	SCS1219--	SCD1219--	12 – 19 cm
Safe-Cuff Small Adult	SCS1725--	SCD1725--	17 – 25 cm
Safe-Cuff Adult	SCS2333--	SCD2333--	23 – 33 cm
Safe-Cuff Adult Long	SCS2333L--	SCD2333L--	23 – 33 cm
Safe-Cuff Large Adult	SCS3140--	SCD3140--	31 – 40 cm
Safe-Cuff Large Adult Long	SCS3140L--	SCD3140L--	31 – 40 cm
Safe-Cuff Thigh	SCS3850--	SCD3850--	38 – 50 cm
Safe-Cuff Vinyl Infant	VSCS0813--	VSCD0813--	08 – 13 cm
Safe-Cuff Vinyl Child	VSCS1219--	VSCD1219--	12 – 19 cm
Safe-Cuff Vinyl Small Adult	VSCS1725--	VSCD1725--	17 – 25 cm
Safe-Cuff Vinyl Adult	VSCS2333--	VSCD2333--	23 – 33 cm
Safe-Cuff Vinyl Adult Long	VSCS2333L--	VSCD2333L--	23 – 33 cm
Safe-Cuff Vinyl Large Adult	VSCS3140--	VSCD3140--	31 – 40 cm
Safe-Cuff Vinyl Large Adult Long	VSCS3140L--	VSCD3140L--	31 – 40 cm
Safe-Cuff Vinyl Thigh	VSCS3850--	VSCD3850--	38 – 50 cm

**CAS Medical Systems, Inc.**  
 44 East Industrial Road  
 Branford, CT 06405 USA  
 Telephone: 203 488 6056  
 Toll Free: 800 227 4414  
 www.casmed.com

--To order cuffs for a specific monitor, add the appropriate suffix to the product number.

For CASMED® 740/750, add **ML** or **HP**

For Critikon Dinamap®, add **DM**

For Philips®, add **HP**

For GE Healthcare®, add **MQ**

For Datascope®, Dräger®, Spacelabs®, add **FL**

For Datex-Ohmeda®, add **MQMF**

Dealer Imprint



ISO 13485  
 Certificate Number FM 72259